

Invacare® Pronto® M91™ Power Wheelchair



Yes, you can.®

Pronto M91™



Anything else is a compromise.™

The Invacare® Pronto® M91™ power wheelchair keeps life simple by offering the flexibility to meet mild to moderate positioning needs all in a base that tackles the more active daily environments. Built on the foundation of SureStep® suspension and true center-wheel drive technology, the Pronto M91 offers proven, durable construction and best-in-class technologies that provide exceptional value.

Center-Wheel Drive



Center wheel drive offers the tightest turning radius possible, as well as improved tracking and ease of driving over all other configurations.

SureStep®



This patented automatic system ensures that everyday obstacles and thresholds are handled -elegantly, for a smooth ride.

Stability Lock



The most pervasive issue with more maneuverable chairs is stability, but this patented technology helps ensure that the chair remains stable, with the consumer in control.

TCD™ Traction Control Design



TCD improves traction in forward and reverse -especially important for transitioning onto ramps and backing up in grass.

MK5™ Electronics



MK5 electronics, the leading power wheelchair electronics developed and designed by Invacare, includes a NX™-75 (non-expandable controller) on the van and rehab seat models with SPJ™+ (Single Purpose Joystick Plus). This combination allows for full turning capability at slow speeds and provides push button on/off, speed selection and horn as well as a bright LED display.

4-Pole Motors



Proven Invacare® 4-pole motors provide high torque and true gear reduction for many rehab needs.

Invacare® Pronto® M91™



Pronto M91 with Van Seat

- Semi-reclining back (standard) for additional comfort
- Choice of cloth or nylon upholstery
- Width adjustable flip-back arms (standard) with optional height adjustment
- Ergonomic seating contours provide support for longer durations of use
- Adjustable height headrest for head and neck support
- Built-in seat pouch on back



Pronto M91-R with Adjustable Rehab Seat


- Adjustable seat and back angle (ASBA) rehab seat is width and depth adjustable in 1" increments.
- Width adjustment range: 16"-20" and 20"-24"
- Depth adjustment range: 16"-19" and 19"-22"
- One seat pan size within each adjustment range
- Angle adjustable back offers range 85° to 105°
- Sling back upholstery provides comfort and support for mild positioning needs
- Round canes accept other back options for additional positioning



Pronto M91 with PTO Plus (Power Tilt Only)

- Adjustable seat and back angle (ASBA) rehab seat is width and depth adjustable in 1" increments
- 18" seat to floor height for easy access to everyday environments
- Up to 55° of tilt for pressure relief (0° - 55°, 5° - 60°)
- Standard with MK5™ NX™-80 controller and SPJ™ 80 (single purpose joystick)
- Weight capacity of 275 lb. or optional 375 lb. with heavy-duty motor package
- Easy battery access through a simple flip back design

Van Seat Upholstery




Grey cloth




Grey vinyl




Tan vinyl



Candy Red



Deep Blue



Magic Black

Shroud Colors


Yes, you can.®
Specifications**

Overall dimensions	Van seat	Adjustable ASBA Rehab Seat	PTO+ Tilt Only
Lengths (with footboard up)	39"		
Lengths (with footboard down)	40" - 42"		
Width (excluding joystick)	25.9"		
Heights	37" - 48"		39"
Seat dimensions			
Widths	18" - 22"	16" - 24"	16" - 22"
Depths	16" - 20"	16" - 22"	16" - 22"
Overall height	37" - 40"	35" - 43"	
Back height	18"	16" - 20"	16" - 24"
with headrest	25"	n/a	n/a
Back angle	80° - 115°	85° - 105°	85° - 105°
Seat-to-floor height	21" - 23"	18.5" - 21"	18" (at 0° of tilt)
Weight			
With batteries	273 lb.		296 lb.
Without batteries	199 lb.		222 lb.
Wheel and caster sizes			
Front casters	6" x 2"		
Drive wheels	14" x 3"		
Rear casters	6" x 2"		
Maximum speed	6.5 mph / 4.5 mph Heavy-duty		
Incline capability	9°		
Weight capacity	300 lb. (400 lb option)		275 lb. (375 lb. option)
Turning radius			
Front (with footboard down)	19.5"		
Rear	21.5"		
Ground clearance	3.25"		
Limited warranties			
Base frame	60 months		
Seat frame	60 months		
Electronics	12 months		
Motor/gearbox	18 months		

* Height measured without headrest and with 20° back

** Due to manufacturing variances, all dimensions, weights and degrees may be $\pm .5"$, $\pm 10\%$ or $\pm 3^\circ$

Invacare Corporation
www.invacare.com

USA

One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada

570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

Invacare, the Blue Medallion Design, Yes, you can., and Pronto are registered trademarks and M91, MK5, NX, SPJ and SureStep are trademarks of Invacare Corporation. Specifications are subject to change without notification.

© 2005

Invacare Corporation
Form No. 02-083
Rev. 10/05