



Seat Colors
Black Vinyl

Body Color
Red

GP160

GP162

**BASE
DISASSEMBLES!**

Roll

**The Ultimate
Retail Cash
Power Chairs!**

GOLDEN

LiteRider[®]ptc

The LiteRider PTC personal transport chairs offer users two great cash retail power chair choices! Both models feature a unique tubular frame design, a comfortable, stadium style seat and 18Ah batteries for an operating range up to 9.5 miles. What makes them different? While both models allow for the seat and battery pack to be easily removed, the GP162 takes disassembly a step further with a base that disassembles quickly and easily!



GOLDEN
technologies

YOUR LIFE IN MOTION

Model Number	GP160	GP162
Weight Capacity	300 lbs.	300 lbs.
Drive Wheels	Rear	Rear
Maximum Speed	3.5 mph	3.5 mph
Operating Range ¹	9.5 miles	9.5 miles
Ground Clearance	2.5"	2.5"
Turning Radius	26"	26"
Type Batteries	2-18 AH	2-18 AH
Colors	Red	Red
Freewheel Mode	Yes	Yes
Electronic Speed Control	Yes	Yes
Electro-Mechanical Brakes	Yes	Yes
Charger	Off Board	Off Board
Controller	PG VR-2 50 Amp Remote	PG VR-2 50 Amp Remote
Dimensions:		
Length	33"	33"
Width	22.5"	22.5"
Height (ground to top of back of seat)	33" - 34"	33" - 34"
Height (ground to top of back of headrest)	N/A	N/A
Ground to Top of Seat	20.5" - 21.5"	20.5" - 21.5"
Number of Seat Height Positions	2	2
Size of Increments	1.0"	1.0"
Weight of Unit: (assembled)		
Base	67 lbs.	N/A
Front Half	N/A.	34 lbs.
Rear Half	N/A	35 lbs.
Seat w/Arms	22 lbs.	22 lbs.
Batteries	25 lbs/pack	25 lbs/pack
Tires:		
Drive Wheels	9"	9"
Casters	6"	6"
Anti-Tip Wheels	3"	3"
Standard Seat:		
Slide Seat	N/A	N/A
Back Height - no headrest	16"	16"
Back Height with headrest	N/A	N/A
Width x Depth	17" x 16"	17" x 16"
Color	Black Vinyl	Black Vinyl
Footrest:		
Height adjustable	Yes	Yes
Angle adjustable	Yes	Yes
Warranty:		
Frame	Five Years	Five Years
Drive Train	13 Months	13 Months
Electronics ²	13 Months	13 Months

¹ Battery range will vary due to rider weight, drive surface and terrain.

² Electronics warranty excludes batteries.

Disassembles easily for transport!



1. Remove seat



2. Pull out battery pack



3. GP162 model only. Grasp the frame and locking handle; lift the handle to unlock, and lift the frame up and off the drive train.



Perfect for transport!

Model GP162 shown.

Available Accessories:

- Cane Holder Tube
- Oxygen Tank Holder
- Cup Holder
- Quad Cane Holder
- Walker Holder
- Crutch Holder
- Safety Flag
- Power Chair Cover
- Touch-up Paint
- Deluxe Pack N' Go
- Seat Belt (Optional)

Some accessories require a mounting bracket or mounting clips.

Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.

