

HR-FLD-115

Multi-Function Deluxe Headrest



- Headrest component can rotate 360, up/down & back or forth.
- Allows compact folding while still attached to the wheelchair.
- Can fit from 16"-20" back width on 7/8" diameter frame.
- Lightweight frame & sturdy padding.
- Easy to maintain.
- Reaches heights from seat 18" to 23"

How To Use

When using wheelchair, unfold the wheelchair first, then prop up the headrest. Conversely when folding the wheelchair, first pull down the headrest to release the support tension, then fold wheelchair.



KARMAN HEALTHCARE Inc.

19255 San Jose Avenue . City of Industry . CA 91748-1418

Order.1-800-805-2762 Tel.626-581-2235 Fax.626-581-2335