

INVACARE®

# Pronto® M41™

Power Wheelchair



Yes, you can.™

The **Invacare® Pronto® M41™** power wheelchair is the latest offering in the sleek and versatile Pronto series of power wheelchairs. The M41 wheelchair features the latest technology in electronics and in-line motors. This combination provides increased efficiency and battery range, as well as an increased top speed. The M41 wheelchair also features a center-wheel drive platform, simplified Invacare® SureStep® Suspension and six wheels on the ground providing stability and maneuverability similar to other Pronto models. All of these features are packed into a simple, clean looking and easy-to-service compact frame that provides mobility with style for people on the go.

### Features

- Compact size for a tight turning radius
- High efficiency in-line motors for greater speed and longer battery range
- SureStep suspension with six wheels on the ground for smooth driving and optimum stability
- Charger connection conveniently located on joystick
- Stylish design
- 3-amp off-board charger
- Three seat options
  - ✓ Medium-back fold-down seat
  - ✓ High-back semi-recline office-style seat
  - ✓ High-back semi-recline van seat with solid seat pan
- Auto style positioning strap
- U1 batteries
- 300 lb. weight capacity

### Pronto® M41™ Wheelchair Specifications

<b>Overall dimensions</b> (measured with fold-down seat)			
Length	34"		
Width (including joystick)	24"		
Height	33"		
Seat-to-floor	19"-23"		
<b>Seat dimensions</b>			
	(fold-down seat)	(semi-recline seat)	(semi-recline seat w/solid seat pan)
Width	18"	16" - 20"	16" - 20"
Depth	17"	16" - 20"	16" - 20"
Seat angle	92°	80° - 115°	90° - 115°
<b>Weight capacity</b>	300 lb.		
<b>Wheel sizes</b>			
Front/rear casters	6" x 2"		
Drive wheels (Flat-free)	10.5" x 3.5"		
<b>Maximum speed</b>	5 mph		
<b>Motor</b>	250 watts		
<b>Incline capability</b>	7° grade		
<b>Turning radius</b>			
Front	17.75"		
Rear	19.5"		
<b>Weight</b>			
	(fold-down seat)	(semi-recline seat)	(semi-recline seat w/solid seat pan)
with U1 batteries	161 lb.	175 lb.	180 lb.
w/out batteries and foot plate	110 lb.	124 lb.	129 lb.
<b>Ground clearance</b>	2"		
<b>Battery</b>	U1		
<b>Battery charger</b>	3 amp off-board		
<b>Per-charge battery range</b>	Up to 17 miles*		
<b>Limited warranties</b>			
Base frame	60 months		
Seat frame	36 months		
Electronics	12 months		
Motor/gearbox	12 months		

\* Travel range will vary depending on user weight, speed, battery charge and terrain. Do not operate on roads, streets or highways.



### Colors Available

Candy Red



Deep Blue



*Pronto M41 wheelchair with Fold-Down Seat*



*Pronto M41 wheelchair with Semi-Recline Office-Style Seat*



*Pronto M41 wheelchair with Solid Seat Pan*



**Invacare Corporation**  
www.invacare.com

### USA

One Invacare Way  
Elyria, Ohio  
44035-4190  
(800) 333-6900

### Canada

570 Matheson Blvd. E.,  
Unit 8  
Mississauga, Ontario  
L4Z 4G4 Canada  
(800) 668-5324

©2010

Invacare Corporation. All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications are subject to change without notification. Form No. 06-033 Rev. 08/10 100633

